

EXPERIMENT NO. 10

Aim : To find the relation between the resonance length and tension of a wire for a constant frequency, using a sonometer.

Apparatus : Sonometer, tuning fork, weights, rubber cork.

Theory : The resonant frequency of a wire of length l is given by :

$$v = \frac{1}{2l} \sqrt{\frac{T}{\mu}} ; \text{ where } T \text{ is the tension and } \mu \text{ its mass/length.}$$

For constant frequency and for a given wire, v and μ are constants.

i.e. the ratio $\frac{\sqrt{T}}{l}$ is a constant and therefore the graph between T & l^2 will be a straight line.

Procedure:

1. Make the pulley of the sonometer frictionless and put a suitable mass, say 1kg, on the hanger.
2. Place a paper rider on the wire to see the vibrations of the wire.
3. Strike the tuning fork on the rubber cork and place it vertically on the sonometer.
4. Adjust the length of the wire by moving the wooden bridges so that the wire & the paper rider vibrates the maximum.
5. Measure the length of the wire between the wooden bridges using the scale attached.
6. Repeat the experiment for masses 1.5kg & 2kg and find the corresponding resonance lengths.
7. Plot a graph between T and l^2


Result : From the graph, we conclude that, $T \propto l^2$

Precautions :

1. Pulley should be frictionless
2. Weight of the hanger should be included

Sources of error :

1. Bridges may not be sharp
2. Pulley may have friction.



Sr. No	Mass on the hanger m (kg)	Tension on the string $(T) = mg$	Resonance length of the wire $l = (l_2 - l_1)$ cm	l^2 (cm^2)
1	0.5			
2	1.0			
3	1.5			

Graph :

